



POLY-GA IMMUNOHISTOCHEMISTRY IS A RELIABLE TOOL FOR DETECTING C9ORF72 GENETIC MUTATION AT THE SYDNEY BRAIN BANK

Author (presenter) Jordan Carroll¹, Heather McCann¹, Glenda Halliday², John BJ Kwok², Claire Shepherd^{1,3}

¹ Neuroscience Research Australia, Barker Street, Randwick NSW 2031

² Faculty of Medicine and Health School of Medical Sciences, University of Sydney Brain and Mind Centre, Camperdown NSW 2050

³ School of Medical Sciences, University of New South Wales, Kensington NSW 2052



THE UNIVERSITY OF
SYDNEY



FRONTIER
Centre for Ageing Research
Frontotemporal Dementia and Motor
Neurodegenerative Syndromes



BACKGROUND

The most common genetic mutation in familial frontotemporal lobar degeneration (FTLD) and motor neuron disease (MND) is a hexanucleotide repeat expansion (G_4C_2) in the *C9orf72* gene. In 2013 a characteristic star-shaped neuronal inclusion pathology was found to contain dipeptide repeat proteins, including poly Gly-Ala (polyGA), which was later found to be the dominant form[1]. The commonly-used repeat-primed PCR method of testing has been shown to have limitations in sensitivity, specificity and inter-laboratory interpretation and not all brain banks perform routine genetic screening of their cases[2]

A To determine the accuracy of Poly-GA
I immunohistochemistry when detecting
M dipeptide repeat pathology in the
S *C9orf72* gene.

METHODS

- Formalin Fixed paraffin embedded human cerebellum tissue from FTLD and MND cases (n=63, n=64).
- Using a commercially available poly-GA antibody (Millipore, MABN889, clone 5E9) at a concentration of 1:1000uL we immunohistochemically stained the cases using a Roche Benchmark GX autostainer in a research setting.
- Diaminobenzidine (DAB) was used as the substrate, and counterstained using Haematoxylin.




Figure 1: Star shaped neuronal cytoplasmic Poly-GA inclusions in the cerebellar cortex[3].



RESULTS




Figure 2: Star and ring shaped neuronal cytoplasmic Poly-GA inclusions in the cerebellar cortex.

- We found 100% accuracy confirming cases with known C9orf72 repeat expansion (n=26).
- A further five cases were found to be positive that had not yet undergone clinical genetic screening. Two cases were subsequently screened and demonstrated a positive result.
- Of the last three cases that did not undergo genetic screening, one had a sibling with C9orf72 mutation and the other two had no known genetic testing results or family history.

PolyGA Positive



■ Known +ve ■ Unknown +ve

Figure 3: Representation of previously unknown +ve C9orf72 mutation.




Figure 4: Star and ring shaped neuronal cytoplasmic Poly-GA inclusions in the granular cell layer of the cerebellum.

CONCLUSION

- PolyGA immunohistochemistry of postmortem human brain tissue serves as a useful tool to accurately identify the presence of dipeptide repeat pathology due to defects in the *C9orf72* gene.
- We have demonstrated the polyGA antibody to be a reliable and sensitive method of detection in a research setting, with high utility in cases where there has been no opportunity for genetic testing during the life of the patient.

REFERENCES

1. Mackenzie, I.R., et al., *Dipeptide repeat protein pathology in C9ORF72 mutation cases: clinico-pathological correlations*. Acta Neuropathol, 2013. **126**(6): p. 859-79.
2. Crook, A., et al., *The C9orf72 hexanucleotide repeat expansion presents a challenge for testing laboratories and genetic counseling*. Amyotroph Lateral Scler Frontotemporal Degener, 2019. **20**(5-6): p. 310-316.
3. Mackenzie, I.R., P. Frick, and M. Neumann, *The neuropathology associated with repeat expansions in the C9ORF72 gene*. Acta Neuropathol, 2014. **127**(3): p. 347-57.



@_GraceWei



grace.wei@sydney.edu.au